

## MO-33-2LFA



Motor sensorless, with Double light (LED), with fluid, sterilizable.  
 Application on General Dentistry and Podiatry.



### Motor data

Dassym reference	106.15.000
Speed Range	1000 up to 40000 rpm
Torque at maximum allowed current	3.5 Ncm

### Dimensions

External diameter	21.8 mm
Motor length	44.7 mm
Weight	81 g

### Other specification

Number of pole pairs	1
Number of phases	3
Type	Brushless, internal spray
Warranty	3 years
Anti-return O-ring	Yes
Sterilization	135 °C
Number of wire connected to the motor	5