

MO-33-LFA



Motor sensorless, with light (LED), with fluid, sterilizable.
 Application on General Dentistry and Podiatry.



Motor data

Dassym reference	106.02.000
Speed Range	1000 up to 40000 rpm
Torque at maximum allowed current	3.5 Ncm

Dimensions

External diameter	21.8 mm
Motor length	44.7 mm
Weight	81 g

Other specification

Number of pole pairs	1
Number of phases	3
Type	Brushless, internal spray
Warranty	3 years
Anti-return O-ring	Yes
Sterilization	135 °C
Number of wire connected to the motor	5