

MO-93-2LFA



Motor with Hall sensors, with Double lights (white+blue/UV), with fluid, sterilizable. Application on General Dentistry and Endodontics.



Motor data

Dassym reference	107.03.000
Speed Range	80 up to 40000 rpm
Torque at maximum allowed current	3.5 Ncm

Dimensions

External diameter	21.8 mm
Motor length	44.7 mm
Weight	81 g

Other specification

Number of pole pairs	1
Number of phases	3
Type	Brushless, internal spray
Warranty	3 years
Anti-return O-ring	Yes
Sterilization	135 °C
Number of wire connected to the motor	5